

8.3 VS 8.2

Clarity (Lite)

ENG

Code/Rev.: M235/80B
Date: 4/14/2020

Phone: +420 251 013 400
Fax: +420 251 013 401
clarity@dataapex.com
www.dataapex.com

DataApex Ltd.
Petrzilkova 2583/13
158 00 Prague 5
The Czech Republic

Sections of the manual connected only to the **Clarity Full** version are marked with the  icon.

Clarity[®], DataApex[®] and ▲[®] are trademarks of DataApex Ltd. Microsoft[®] and Windows[™] are trademarks of Microsoft Corporation.

DataApex reserves the right to make changes to manuals without prior notice. Updated manuals can be downloaded from www.dataapex.com.

Author: MP

Contents

1 Preamble	1
2 Clarity	2
2.1 Method versioning	2
2.1.1 General concepts	2
2.1.2 Audit Trail	3
2.1.3 Extended functionality of Save as...	3
2.2 Calibration	5
2.2.1 Grouped parameters on the tab of a specific compound	5
2.2.2 Removed Weight from Calibration Options	6
2.2.3 Improved visibility of recalibrations	6
2.3 Audit Trail	7
2.4 Extension MS	8
2.4.1 MS Method	8
2.4.2 MS Integration	8
2.4.3 Report Setup	10
2.5 Other changes	11
3 New and updated control modules	12
3.1 Agilent	12
3.2 Analytik Jena	12
3.3 Antec	12
3.4 CQS	12
3.5 Dani	12
3.6 Esensing	13
3.7 Fuli	13
3.8 JAI	13
3.9 Futecs	13
3.10 Knauer	13
3.11 Sedere	13
3.12 Shimadzu	13
3.13 Spark	13
3.14 Sykam	14
3.15 VICI Valco Instruments	14
3.16 Watrex	14
3.17 Young In Chromass	14

To facilitate the orientation in the **8.3 vs 8.2** manual and **Clarity** chromatography station, different fonts are used throughout the manual. Meanings of these fonts are:

Instrument (blue text) marks the name of the window to which the text refers.

Open File (italics) describes the commands and names of fields in **Clarity**, parameters that can be entered into them or a window or dialog name (when you already are in the topic describing the window).

WORK1 (capitals) indicates the name of the file and/or directory.

ACTIVE (capital italics) marks the state of the station or its part.

The bold text is sometimes also used for important parts of the text and the name of the **Clarity** station. Moreover, some sections are written in format other than normal text. These sections are formatted as follows:

Note: Notifies the reader of relevant information.

Caution: Warns the user of possibly dangerous or very important information.

Marks the problem statement or trouble question.

Description: Presents more detailed information on the problem, describes its causes, etc.

Solution: Marks the response to the question, presents a procedure how to remove it.

1 Preamble

This document will guide you through the news and improvements in the **Clarity** Chromatography Station version **8.3** compared to version **8.2**.

The most interesting features of version 8.3 include:

- Method versioning
- Changes in the Calibration - grouped parameters on the tab of a specific compound
- Audit Trail - new columns and icons
- MS Extension enhancements
- New and updated control modules

2 Clarity

2.1 Method versioning

2.1.1 General concepts

It is sometimes necessary to distinguish with what method parameters was a chromatogram acquired. Therefore we have implemented method versioning, just like in the chromatogram.

Method versioning is intended for tracking down changes in the method file (*.met). With each save, a new version of the method is created. These method versions are saved internally into the corresponding method file.

Any existing method (from any previous version of Clarity) will now display version #1, even though it may have been saved with new parameters numerous times before, and now any subsequent method save will create a new version - incrementing the version by +1.

The power of method versioning is that you can go back in time and see what were the parameters at any specific time.

[Select Method](#) dialog invoked from the [Method Setup](#) window (*Open...* command) now contains *Version* listbox allowing to open any historical version of the method file.




Fig 1: Select Method

To easily track what method version is opened, [Method Setup](#) header now displays version number and time of last method save.




Fig 2: Method Setup

Note: Any historical version of the method will be opened in read-only mode which means it is not possible to change any parameters. To edit historical version you must save it under a new name.

2.1.2 Audit Trail

Changes to the method are logged into the **Method Audit Trail**, along with the information about the version and save date.

Time	Group	Severity	User Name	Det	Description
03.12.2019 13:51:21	File		Administrator		Save File C:\ClarityREQ2016\DataFiles\DEMO1\Demo1.met - #2; 03.12.2019 13:51:21
03.12.2019 13:51:21	Method		Administrator		Measurement: Run Time has been changed from 1,7 to 2
03.12.2019 13:50:39	File		Administrator		Save File C:\ClarityREQ2016\DataFiles\DEMO1\Demo1.met - #1; 03.12.2019 13:50:39
03.12.2019 13:50:38	Method		Administrator		DEMO - 4: Range was changed from 12000 to 12500 mV.
03.12.2019 13:50:38	Method		Administrator		DEMO - 3: Range was changed from 12000 to 1250 mV.
03.12.2019 13:50:38	Method		Administrator		DEMO - 2: Range was changed from 10000 to 12000 mV.
03.12.2019 13:50:38	Method		Administrator		DEMO - 1: Range was changed from 10000 to 12000 mV.
03.12.2019 13:50:37	Method		Administrator		Method adapted - This method was adapted to suit new instrument configuration.
14.06.2018 15:00:00			Administrator		Save File Noname.met As C:\Clarity\DataFiles\DEMO1\Demo1.met

Fig 3: Method Setup

2.1.3 Extended functionality of Save as...

Along with method versioning, the **Save As...** functionality has been extended. The default remains the same, clicking directly the **Save As...**

icon invokes dialog to specify the new file name, after that the newly created method is loaded in the [Method Setup](#) dialog.

Clicking the small arrow next to the *Save As...* icon reveals the available options:

- Makes a copy of the method along with any changes and opens it (changes in the original method are discarded).
- Makes a copy of the method along with any changes but does not open it and you continue working with the original one.




Fig 4: Method Setup

2.2 Calibration

2.2.1 Grouped parameters on the tab of a specific compound

From the identification and quantification point of view, **Calibration** is perhaps the most important part of data processing. We have redesigned the tab of a specific compound to better distinguish individual dependencies and grouped them into logical units.




Fig 5: Calibration - tab of a specific compound

The left part of **Calibration** window is now divided into 3 groups (table with levels is ignored):

- ① Identification - it groups parameters which define conditions under which the peak is recognized.

- ② Quantification - it groups parameters which affect the calculation of the calibration curve.
- ③ Calibration curve - it groups parameters which are calculated from the calibration curve. All parameters displayed here are read-only and cannot be changed.

Some parameters in the groups above are context based, meaning changing the value in one field directly affects other field. Such example is the *Curve Fit Type* which affects whether the option below it will be *Origin*, *Manual Response Factor* or *Ending Point*.

Some Clarity Extensions feature additional amendments but the principle remains the same.

2.2.2 Removed Weight from Calibration Options

Weight option was a historical remain which users found confusing and it was thus removed from the [Calibration Options](#). Any calibration set to *Weight* option will be automatically transferred to *Average* and the change will be recorded in the *Calibration Audit Trail*.

The screenshot shows the 'Calibration Options (Calib\250X8HR1)' dialog box. It has a 'Defaults' tab selected. The 'Calibration Description' field contains 'IEX H Form, 9 mM H2SO4, 0.5 ml/min'. The 'Display Mode' is set to 'ESTD'. The 'Number of Signals' is 1. Under 'Calibration', 'Automatic' is selected. Under 'Apply on', 'On All Signals' is selected. 'Compound Units' is 'g/l'. Under 'Recalibration', 'Replace' is selected and 'No. of Points' is 10. There are checkboxes for 'Recalibration Search Criteria' (unchecked, 1%), 'Enable Manual Response Value Change' (checked), 'Update Retention Time' (checked), 'Default Injected Volume' (unchecked, 0 µl), 'Retention Indexes use Log. Interpolation with Unretained Peak' (checked), and 'Response Factor as Response / Amount' (unchecked). The 'Curve Check' section has 'Deviation' (0%) and 'Correlation' (0%) fields. At the bottom are 'OK', 'Cancel', and 'Help' buttons.

Fig 6: Calibration Options

2.2.3 Improved visibility of recalibrations

Rec No. ④ (Fig 5 on pg 5.) shows number of used recalibrations / total recalibrations on the respective level. Clicking ... invokes the [Details of Calibration Point](#) dialog with a detailed log of the selected level of the calibration point. Levels that have not yet been (re)calibrated have default value 0/0.

2.3 Audit Trail ✓ Full version

Audit Trail serves as a log of individual operations of the station. It is required by the GLP practice. It is also useful for onsite troubleshooting and technical support where we can see what was happening on the station.

We have added two more columns:

- **Group** column displays icon and group into which the given audit trail entry falls.
- **Severity** column displays the severity of the logged operation. Severity can be either <empty>, *Warning* or *Error*.

In the future we plan to allow filtering that will enable to easily search for vital information. Therefore the "cleanup" was intended to make better sense of the logged operations.

OK	Time	Group	Severity	User Name	Instrument	Area	Description	Info
82	06.11.2019 08:59:07	File		Administrator	Instrument	Files	Open File C:\Clarity833\DataFiles\DEMO1\CalmExample.caf	Clarity 8.3.0.62 FULL; SN: 088-000888
83	06.11.2019 08:59:07	File		Administrator	Instrument	Acquisition	Running - Acquisition running; Acquisition was stopped by detection of stop of all detectors	Clarity 8.3.0.62 FULL; SN: 088-000888
84	06.11.2019 08:59:07	Stop		Administrator	Instrument	Acquisition	Acquisition was stopped by detection of stop of all detectors	Clarity 8.3.0.62 FULL; SN: 088-000888
85	06.11.2019 08:57:25	Start		Administrator	Instrument	Acquisition	Waiting - Waiting for injector; Acquisition was started by digital input	Clarity 8.3.0.62 FULL; SN: 088-000888
86	06.11.2019 08:57:24	Injection		Administrator	Instrument	Acquisition	Waiting - Waiting for injector; Perform Injection	Clarity 8.3.0.62 FULL; SN: 088-000888
87	06.11.2019 08:57:23	Acquisition		Administrator	Instrument	Detector	Colbrick - 4 (SN 1284); Method sent to sub-device.	Clarity 8.3.0.62 FULL; SN: 088-000888
88	06.11.2019 08:57:23	Acquisition		Administrator	Instrument	Detector	Colbrick - 3 (SN 1284); Method sent to sub-device.	Clarity 8.3.0.62 FULL; SN: 088-000888
89	06.11.2019 08:57:23	Acquisition		Administrator	Instrument	Detector	Colbrick - 2 (SN 1284); Method sent to sub-device.	Clarity 8.3.0.62 FULL; SN: 088-000888
90	06.11.2019 08:57:23	Acquisition		Administrator	Instrument	Detector	Colbrick - 1 (SN 1284); Method sent to sub-device.	Clarity 8.3.0.62 FULL; SN: 088-000888
91	06.11.2019 08:57:23	Acquisition		Administrator	Instrument	Acquisition	Sending method - Sending method to hardware; Method was sent into instrument	Clarity 8.3.0.62 FULL; SN: 088-000888
92	06.11.2019 08:57:23	Acquisition		Administrator	Instrument	Acquisition	Method C:\Clarity833\DataFiles\DEMO1\Demo1.met Ver. 1.06.11.2019 08:57:17 (Size: 325063, CRC: 979)	Clarity 8.3.0.62 FULL; SN: 088-000888
93	06.11.2019 08:57:23	File		Administrator	Instrument	Files	Open File C:\Clarity833\DataFiles\DEMO1\CalmExample.caf	Clarity 8.3.0.62 FULL; SN: 088-000888
94	06.11.2019 08:57:23	Acquisition		Administrator	Instrument	Sequence State	Injecting Row=1, Val=# 1 (1), Trj=1, SampleID=146030000, Volume=5,000[µL], File=#02, Sequence=	Clarity 8.3.0.62 FULL; SN: 088-000888
95	06.11.2019 08:57:23	Sequence		Administrator	Instrument	Sequence State	Sequence C:\Clarity833\DataFiles\DEMO1\Demo1.met started C:\Clarity833\DataFiles\DEMO1\Demo1.met	Clarity 8.3.0.62 FULL; SN: 088-000888
96	06.11.2019 08:57:23	Sequence		Administrator	Instrument	Acquisition	Ready - Ready to start run; Sequence was started by user	Clarity 8.3.0.62 FULL; SN: 088-000888
97	06.11.2019 08:57:23	Sequence		Administrator	Instrument	Sequence	C:\Clarity833\DataFiles\DEMO1\Demo1.met - Sequence Start Conditions: Current Value of Counter 1.	Clarity 8.3.0.62 FULL; SN: 088-000888
98	06.11.2019 08:57:23	Sequence		Administrator	Instrument	Sequence	C:\Clarity833\DataFiles\DEMO1\Demo1.met - Sequence Start Conditions: Lines to Run 1-5.	Clarity 8.3.0.62 FULL; SN: 088-000888
99	06.11.2019 08:57:23	File		Administrator	Instrument	Files	Open File C:\Clarity833\DataFiles\DEMO1\CalmExample.caf	Clarity 8.3.0.62 FULL; SN: 088-000888
100	06.11.2019 08:57:23	File		Administrator	Instrument	Files	Open File C:\Clarity833\DataFiles\DEMO1\CalmExample.caf	Clarity 8.3.0.62 FULL; SN: 088-000888
101	06.11.2019 08:57:19	File		Administrator	Instrument	Files	Open File C:\Clarity833\DataFiles\DEMO1\CalmExample.caf	Clarity 8.3.0.62 FULL; SN: 088-000888
102	06.11.2019 08:57:19	File		Administrator	Instrument	Files	Open File C:\Clarity833\DataFiles\DEMO1\CalmExample.caf	Clarity 8.3.0.62 FULL; SN: 088-000888
103	06.11.2019 08:57:19	File		Administrator	Instrument	Files	Open File C:\Clarity833\DataFiles\Common\Sequence.sty	Clarity 8.3.0.62 FULL; SN: 088-000888
104	06.11.2019 08:57:18	File		Administrator	Instrument	Files	Open File C:\Clarity833\DataFiles\DEMO1\Demo1.met	Clarity 8.3.0.62 FULL; SN: 088-000888
105	06.11.2019 08:57:17	Acquisition		Administrator	Instrument	Detector	Colbrick - 4 (SN 1284); Method sent to sub-device.	Clarity 8.3.0.62 FULL; SN: 088-000888
106	06.11.2019 08:57:17	Acquisition		Administrator	Instrument	Detector	Colbrick - 3 (SN 1284); Method sent to sub-device.	Clarity 8.3.0.62 FULL; SN: 088-000888
107	06.11.2019 08:57:17	Acquisition		Administrator	Instrument	Detector	Colbrick - 2 (SN 1284); Method sent to sub-device.	Clarity 8.3.0.62 FULL; SN: 088-000888
108	06.11.2019 08:57:17	Acquisition		Administrator	Instrument	Detector	Colbrick - 1 (SN 1284); Method sent to sub-device.	Clarity 8.3.0.62 FULL; SN: 088-000888
109	06.11.2019 08:57:17	Acquisition		Administrator	Instrument	Acquisition	Method C:\Clarity833\DataFiles\DEMO1\Demo1.met Ver. 1.06.11.2019 08:57:17 (Size: 325063, CRC: 979)	Clarity 8.3.0.62 FULL; SN: 088-000888
110	06.11.2019 08:57:17	File		Administrator	Instrument	Files	Save File C:\Clarity833\DataFiles\DEMO1\Demo1.met Ver. 1.06.11.2019 08:57:17	Clarity 8.3.0.62 FULL; SN: 088-000888
111	06.11.2019 08:57:17	Method		Administrator	Instrument	Method Change	C:\Clarity833\DataFiles\DEMO1\Demo1.met - Colbrick - 4 (SN 1284); Rate was changed from 10,00 to 12,	Clarity 8.3.0.62 FULL; SN: 088-000888
112	06.11.2019 08:57:17	Method		Administrator	Instrument	Method Change	C:\Clarity833\DataFiles\DEMO1\Demo1.met - Colbrick - 3 (SN 1284); Rate was changed from 10,00 to 12,	Clarity 8.3.0.62 FULL; SN: 088-000888
113	06.11.2019 08:57:17	Method		Administrator	Instrument	Method Change	C:\Clarity833\DataFiles\DEMO1\Demo1.met - Colbrick - 2 (SN 1284); Rate was changed from 20,00 to 25,	Clarity 8.3.0.62 FULL; SN: 088-000888
114	06.11.2019 08:57:17	Method		Administrator	Instrument	Method Change	C:\Clarity833\DataFiles\DEMO1\Demo1.met - Colbrick - 2 (SN 1284); Range was changed from 10000 to 1,	Clarity 8.3.0.62 FULL; SN: 088-000888
115	06.11.2019 08:57:17	Method		Administrator	Instrument	Method Change	C:\Clarity833\DataFiles\DEMO1\Demo1.met - Colbrick - 1 (SN 1284); Rate was changed from 20,00 to 25,	Clarity 8.3.0.62 FULL; SN: 088-000888
116	06.11.2019 08:57:17	Method		Administrator	Instrument	Method Change	C:\Clarity833\DataFiles\DEMO1\Demo1.met - Colbrick - 1 (SN 1284); Range was changed from 10000 to 1,	Clarity 8.3.0.62 FULL; SN: 088-000888
117	06.11.2019 08:57:16	File		Administrator	Instrument	Files	Open File C:\Clarity833\DataFiles\DEMO1\CalmExample.caf	Clarity 8.3.0.62 FULL; SN: 088-000888
118	06.11.2019 08:57:16	Method		Administrator	Instrument	Method Change	C:\Clarity833\DataFiles\DEMO1\Demo1.met - Method adapted - This method was adapted to suit new inst	Clarity 8.3.0.62 FULL; SN: 088-000888

Fig 7: Audit Trail

2.4 Extension MS Full version

We continually strive to improve the user experience. This was the main reason for the changes which are described next.

2.4.1 MS Method

MS Method is for the MS extension an important part as it holds information about the evaluation of spectral data. *MS Method* is no more hidden under the *Measurement Conditions* tab. It is on the same level as *Results*, *Summary* and similar tabs. We hope that this change will enable users to access settings of the MS Method more easily and faster.




Fig 8: Chromatogram - MS Method tab

2.4.2 MS Integration

Integration table is quite context based - information displayed is based according to which signal is selected. However MS Extension adds one more level of confusion - the quantification signal which can also be integrated.

Therefore the *MS Integration* tab has been isolated from the *Integration* tab. In the *MS Integration* tab only quantification signals are displayed and

can be integrated. To change integration of a quantification signal - from the *Quant. Signal* combobox choose which signal you want to integrate and then continue just as you would in the standard *Integration* table.




Fig 9: Chromatogram - MS Integration tab

Integration tab (solely in the MS Extension) now displays only signals which come from the *Result Table*. *Signal* combobox lists detector signals only - without quantification signals.




Fig 10: Chromatogram - Integration tab

2.4.3 Report Setup

MS tab in the Report Setup now contains section MS Integration Table which allows to select which integration table of quantification signal (option to print all signals is available) will be printed.




Fig 11: Report Setup - MS tab

2.5 Other changes

- Audit Trail: AreaEN and DescriptionEN columns will be filled even when the station is switched into English.
- Audit Trail: Audit Trails created with version 8.3 cannot be opened in previous versions due changes in the internal format.
- Chromatogram: Enabled a limited support for Paste operation into the Integration table.
- Chromatogram: Improved Asymmetry calculation in case of too few datapoints.
- GLP: New item "Disallow Chromatogram Merge Operation" in the GLP Options dialog.
- IQ: Reorganized IQ to better distinguish between control modules installed by Clarity and by 3rd parties.
- Installation: Installer step for selecting acquisition device was removed from the installation as it confused many users. This removal will not affect any existing installations.
- Installation: Entering User Code (U/C) step was moved almost to the end. Once the U/C is entered, the installation jumps to the last step where the user may start Clarity or simply end the installation process.
- OQ Validation: Instead of Knauer Virtual Detector a new OQ Validation detector by Ruby script is used.
- PDA Extension: Changes in Peak Purity calculations. The reference spectrum now is peak apex in the spectrum (previously in the active signal, which may have been shifted). In case the calculated correlation is negative, it is replaced by zero (previously absolute value was used). Corrected calculation of absorbance threshold.

3 New and updated control modules

This section contains new and updated control modules introduced in Clarity 8.3.

Testing state is dedicated for new control modules.

Ready state is dedicated for existing control modules that have been updated or somehow improved.

3.1 Agilent

New:

- Agilent 8860, 8890 and Intuvo 900 GC control modules are now in the Testing state.

Updated:

- Agilent ICF libraries updated to version A.02.06.DU3.

Removed:

- Agilent 1200 control module supporting G1315C, G1315D, G1321B detectors was removed from installation. In order to continue using Agilent 1200, please use the Agilent ICF LC driver.

3.2 Analytik Jena

New:

- Analytik Jena PQ LC / IC control module is now in the Testing state.

3.3 Antec

Updated:

- Antec Decade Elite control module - added support for pulse mode 2. See respective manual for further details.
- Antec AS 110 control module was replaced by Alias II. control module - the functionality remained the same.

3.4 CQS

Updated:

- Climax HPLC system control module is now in the Ready state.

3.5 Dani

Updated:

- Dani Master GC control module updated to version 1.6.11.0.

3.6 Esensing

Updated:

- MassChrom1299 HPLC system control module is now in the Ready state.

3.7 Fuli

Updated:

- GC 9720 Plus control module is now in the Ready state.
- GC 9790 Plus control module is now in the Ready state.

3.8 JAI

Updated:

- JAI FC-7000 control module is now in the Ready state.

3.9 Futecs

New:

- Futecs CD-6000, CT-6000, P-6000 and RI-6000 control modules are now in the Testing state.

3.10 Knauer

Updated:

- Knauer HPLC control module updated to version 8.1.0.6036.

3.11 Sedere

Updated:

- Sedere Sedex 85/90 and ELSD USB drivers updated to version 1.3.
- Sedere Sedex LC/FP 100 USB drivers updated to version 1.2.

3.12 Shimadzu

Updated:

- Shimadzu GC2010 control module is now in the Ready state.

3.13 Spark

New:

- Spark LC1299 and LC1299 Plus HPLC systems control modules are now in the Testing state.

Updated:

- Spark Symbiosis ACE - added support for new options for addressing cartridges.

3.14 Sykam

New:

- Sykam S4120 control module is now in the Testing state.
- Sykam S150 control module is now in the Testing state.

Updated:

- Sykam S6510 control module is now in the Ready state.
- Sykam S1130 control module updated to version 2.0.1.5.
- Sykam S5300 control module updated to version 2.0.0.13.

3.15 VICI Valco Instruments

Updated:

- VICI Valco TCD3 Thermal Conductivity detector control module is now in the Ready state.
- VICI Valco valves - option to set position while closing Instrument has been moved from [Method Setup](#) to [System Configuration](#) of the respective valve.

3.16 Watrex

New:

- Watrex Streamline CT1, IC1, P1m and P1s control modules are now in the Testing state.

Updated:

- Watrex Streamline AS1, AS2 and UV1 control modules are now in the Ready state.

3.17 Young In Chromass

New:

- Young In Chromass ChroZen GC control module is now in the Testing state.

Updated:

- Young In Chromass YL9130 control module updated to version 4.0.1.12.
- Young In Chromass YL9150 Plus control module updated to version 1.0.0.6.