

## 8.0 VS 7.4

Clarity (Lite)

ENG

Code/Rev.: M226/80B  
Date: 6/18/2018

Phone: +420 251 013 400  
Fax: +420 251 013 401  
clarity@dataapex.com  
www.dataapex.com

DataApex Ltd.  
Petrzilkova 2583/13  
158 00 Prague 5  
The Czech Republic

Sections of the manual connected only to the **Clarity Full** version are marked with the  icon.

Clarity<sup>®</sup>, DataApex<sup>®</sup> and ▲<sup>®</sup> are trademarks of DataApex Ltd. Microsoft<sup>®</sup> and Windows<sup>™</sup> are trademarks of Microsoft Corporation.

*DataApex reserves the right to make changes to manuals without prior notice. Updated manuals can be downloaded from [www.dataapex.com](http://www.dataapex.com).*

Author: MP

# Contents

<b>1 Preamble</b> .....	<b>1</b>
<b>2 Clarity</b> .....	<b>2</b>
2.1 Main window .....	2
2.2 Instrument window .....	3
2.2.1 The window is divided into 3 separate parts: .....	3
2.2.2 Revised Instrument states .....	4
2.2.3 System Configuration - Method Sending Options .....	6
2.3 Method Setup .....	7
2.4 Single Analysis .....	9
2.5 Device Monitor .....	11
2.6 Icons .....	12
2.7 Batch .....	14
2.8 Good Laboratory Practice (GLP) .....	15
2.9 MS Extension - enhancements .....	16
2.9.1 Unified features .....	16
2.9.2 Single Compound Search .....	16
2.9.3 Automatic Compound Search .....	17
2.9.4 Independent size of views .....	18
2.10 Internal Standards (ISTD) .....	20
2.11 Calibration .....	22
2.11.1 Details of Calibration point .....	22
2.12 Chromatogram .....	23
2.12.1 Enhanced error format .....	23
2.12.2 Import Chromatogram .....	23
2.12.3 Mathematical Operations .....	24
2.12.4 New icons for Set Model and Copy from Model .....	25
2.13 Print .....	26
2.13.1 Report Setup - Word Wrap Long Texts .....	26
2.13.2 Report Setup - Force Inverted .....	26
2.13.3 Report Setup - new Lab. Header variable .....	27
2.13.4 Renamed Report Style Instrument.sty .....	27
2.14 Installation .....	28
2.14.1 Windows XP are no longer supported by Clarity 8.0 .....	28
2.14.2 Changes in the Typical installation .....	28
2.15 Other changes .....	29
<b>3 New and updated control modules</b> .....	<b>30</b>
3.1 Advion .....	30
3.2 CoruiTech .....	30
3.3 CTC .....	30
3.4 Dani .....	30
3.5 Ecom .....	30
3.6 HTA .....	30
3.7 KNAUER .....	31
3.8 Shimadzu .....	31

---

3.9 Schambeck .....	31
3.10 YL Instruments .....	31

To facilitate the orientation in the **8.0 vs 7.4** manual and **Clarity** chromatography station, different fonts are used throughout the manual. Meanings of these fonts are:

**Instrument** (blue text) marks the name of the window to which the text refers.

*Open File* (italics) describes the commands and names of fields in **Clarity**, parameters that can be entered into them or a window or dialog name (when you already are in the topic describing the window).

WORK1 (capitals) indicates the name of the file and/or directory.

*ACTIVE* (capital italics) marks the state of the station or its part.

The bold text is sometimes also used for important parts of the text and the name of the **Clarity** station. Moreover, some sections are written in format other than normal text. These sections are formatted as follows:

---

**Note:** Notifies the reader of relevant information.

---

**Caution:** Warns the user of possibly dangerous or very important information.

---

**■ Marks the problem statement or trouble question.**

**Description:** Presents more detailed information on the problem, describes its causes, etc.

**Solution:** Marks the response to the question, presents a procedure how to remove it.

# 1 Preamble

This document will guide you through the news and improvements in the **Clarity Chromatography Station version 8.0** compared to version **7.4**.

**The most interesting features of version 8.0 include:**


- Graphically enhanced user interface for:
  - Main window
  - Instrument window
  - Method Setup dialog
  - Single Analysis dialog
  - Device Monitor window
  - Icons
- Batch dialog - method related changes
- Improvements to GLP compliance - Global Audit Trail
- Improvements in MS Extension
- Increased number of ISTDs
- Calibration - Show details of Calibration point
- Chromatogram:
  - Improvements to error/warning notification in Result Table and All Signals Result Table
  - Option to choose method in the Import and Mathematical Operations dialogs
- Improvements to Print - Report Setup
- Clarity **8.0** can't be installed on Windows XP
- New and updated control modules

# 2 Clarity

## 2.1 Main window

**Main** window has undergone a major graphical change. From the functional point of view, it remains the same. Meaning it still serves for logging into configured Instruments.

However the window is now more understandable and better looking.



*Fig 1: Main window*

- Each Clarity Instrument is distinguished by its color.
  - In the image above, Instrument called "My GC" is blue and therefore all subsequent windows opened from this Instrument will have this color in the header. This analogy was kept from the previous version where the Instrument had different colors of the pipeline.
  - Opened/Locked Instruments will display a thin stripe in the color of the respective Instrument.
- Second Instrument ("My LC") is not red but dark yellow now.
  - For many users it caused confusion in the Audit Trail as operations logged in red (by the color of the Instrument) looked like an error.
- Instrument 4 is Offline meaning there is insufficient number of purchased licenses.
  - Such Instrument is possible to use for data processing such as evaluation of measured chromatograms, preparation of sequence, but cannot be used for data acquisition nor for instrument control.
- Icons on the left vertical toolbar are now bigger for easier accessibility.
- We have removed the option to switch **Main** Clarity window to so-called Clarity Switch (alternate and much smaller Main window).

## 2.2 Instrument window

**Instrument** window has undergone not only graphical change but a lot of functional changes. This section will walk you through the most important ones.




Fig 2: Instrument window description

### 2.2.1 The window is divided into 3 separate parts:

- Upper part with icons opening: [Method Setup](#), [Single Analysis](#), [Sequence](#), [Device Monitor](#), [Data Acquisition](#), [Chromatogram](#) and [Calibration](#).
- Middle part displaying status of the **Instrument** along with running analysis time.
- Lower part displaying a detailed information about running analysis.
- Major differences between Clarity 8.0 and previous versions:
  - The so-called "Instrument Method" or "Method opened on the Instrument" has been abandoned. For a significant number of users it caused confusion which method was used for analysis. Therefore method that will be used for measurement using the Single Analysis was moved to the [Single Analysis](#) dialog where it makes a better sense.
  - We have removed multiple icons leading to different tabs of the [Method Setup](#). Users could not tell the difference between them so now we have one icon for the [Method Setup](#) dialog. It opens on the last closed tab of the [Method Setup](#) dialog. In the menu there is still the option to open a particular method tab.
  - Work with method such as creating a new one or saving was moved directly to the [Method Setup](#) dialog.
    - [Method Setup](#) icon from the **Instrument** window always opens with the last opened method.



- Method that is opened (displayed) in the **Method Setup** dialog must not necessarily be the same as the one that is being measured by. Information table in the **Instrument** window displays *Sent Method* field which indicates according to which method the measurement is or will be carried out.
- In order to measure by single analysis, choose the desired method in the **Single Analysis** dialog, click the *Send Method* button and then the *Run* button.
- *PostRun Settings* that were located in the **Instrument** window but served only for measurement using the **Single Analysis** were moved to the **Single Analysis** dialog.
- **Instrument** window now supports horizontal and vertical re-sizeability.
  - Horizontally it can be increased as much as it allows
  - Vertical: hiding or showing the bottom part (information table) of the **Instrument** window can be achieved by pressing the arrow separating the instrument status and the information table.

## 2.2.2 Revised Instrument states

Instrument states are indispensable as they provide information about the state of the Instrument. Quite often were users having troubles with the meaning of the Instrument states, therefore we have decided to provide a short (quick overview) and long (informative) description. Table bellow sums up the changes:

Instrument state in 7.4.1 and older	Short Instrument state in 8.0	Long Instrument state in 8.0
Communication Error	Error	Error in communication with hardware
Disabled	Disabled	No configured detector
Disabled	Offline	Offline instrument
Control	Control	Performing remaining time program
Sequence Paused	Paused	Sequence paused
Waiting	Ready	Ready to start run
Ready	Ready	Ready for injection
Sending Method	Sending Method	Sending method to hardware
Not ready (Method has not been sent)	No method sent	Ready to send method or start sequence
Not ready to send method	Waiting	Hardware not ready to receive method

Instrument state in 7.4.1 and older	Short Instrument state in 8.0	Long Instrument state in 8.0
Not ready to run acquisition	Waiting	Hardware not ready to run
Waiting for external start	Waiting	Waiting for external start
Waiting for ready	Waiting	Waiting for hardware to be ready to run
Waiting for ready to receive method	Waiting	Waiting for hardware to be ready to receive method
Waiting for injection	Waiting	Waiting for injection
Idle	Waiting (Idle)	Waiting remaining idle time
Running	Running	Acquisition running
Stopping	Running (Finishing)	Waiting for remaining data from hardware
Running (Stopping)	Running (Stopping)	Acquisition running (stop sequence after current acquisition)
Pausing	Running (Pausing)	Acquisition running (pause sequence after current acquisition)

## 2.2.3 System Configuration - Method Sending Options

Changes connected to the method handling have resulted in changes in the [System Configuration - Method Sending Options](#). The option "Prompt for Confirmation" was removed from the dialog and the default option is "Do Not Send Method to Instrument".

There are two options now:

- *Send Method to Instrument* - automatically sends method (that is displayed in the [Single Analysis](#) dialog) to all connected hardware. This option is not recommended for HPLC systems and should be set only by experienced users.
- *Do Not Send Method to Instrument* - does not automatically send method to all connected hardware. Prior to starting acquisition, the method must be send manually.

*Note:* *Method Sending Options* is not accessible while only DataApex's HW (such as INT9, Colibrick, etc.) is configured on the respective instrument. In that case, the method is being sent automatically.




Fig 3: System Configuration - Method Sending Options

## 2.3 Method Setup

For easier understanding, the manipulation with method was moved from the [Instrument](#) window to the [Method Setup](#) dialog.




Fig 4: Method Setup dialog

The title of the [Method Setup](#) dialog displays a method that is currently opened. In the figure above, it's the DEMO1.

Individual tabs of the method have not changed so there is no compatibility issues when migrating from older versions of Clarity.

### Major differences between Clarity 8.0 and previous versions:

- Newly the upper part of the [Method Setup](#) dialog displays a set of icons with which it's possible to create a new method, open existing method, save method, save method as, open report setup, open method audit trail, send method by e-mail or open help.
- If you make any changes, the method becomes Modified, just as you were used to. Upon pressing the **OK** button it will be automatically saved and [Method Setup](#) dialog will close.

- Only a method that is saved can be sent, it's not possible to send a method that is Modified. In such case the user will be notified. Remedy is to *Save* the method and click on *Send Method* once again.
- Pressing the *Send Method* button will result in two actions:
  - Method will be sent to all connected hardware and thus displayed in the information table of the **Instrument** window - just as before.
  - Method will be set as *Method* for **Single Analysis** - a new feature.

## 2.4 Single Analysis

With the redesigned **Instrument** window some parts were moved into **Single Analysis** dialog where it is logical to place them.

Fig 5: Single Analysis dialog


### Major differences between Clarity 8.0 and previous versions:

- **Single Analysis** tabs:
  - Previous main functionality of the **Single Analysis** dialog such as *Sample*, *Sample ID*, *Amount*, *Sample Type*, etc. was moved onto individual tab called *Analysis*.
  - *Post Run Settings* which were located in the **Instrument** window but the actions were applicable only during measurement of **Single Analysis** - for this reason it's more suitable to have the options here.
  - Instead of being in an expandable section, *User Variables* are now also on a separate tab.
- The biggest change is that user needs to select a method that will be used for the **Single Analysis** measurement. However it should be clear which method is being used for the acquisition - method that is displayed in the *Method* field.
  - Method can be either chosen using the three dot button or by directly typing the method name.

- If you open **Method Setup** dialog, open method and then click *Send Method*, the chosen method will be send to all configured hardware and will be set to the *Method* field on the *Analysis* tab.
- If you want to print results, simply choose the *Report Style* that will be used once *Print* is selected on the *Post Run Settings* tab.
- *Sample Type*, *Method* and *Report Style* fields will be disabled during running acquisition meaning that they cannot be changed. To change them you need to either *Stop*, *Abort* or let the analysis to finish.

## 2.5 Device Monitor

**Device Monitor** has changed only graphically to better suit the rest of Clarity.



*Fig 6: Device Monitor window*

**Note:** Even though that we would like to have all the **Device Monitor** windows for all control modules the same, in the sense of unified graphical user interface, with externally developed control modules, it's not possible. It may happen that some sections of the **Device Monitor** window will have older (blue background) design which should not influence the functionality of the control module. This issue will be fixed in the coming versions.



## 2.6 Icons

Clarity 8.0 contains a whole new set of icons which are better looking, more intuitive and responsive (i.e. will increase in size based on the screen resolution).

If you have in any way customized (other than default) your toolbars it is necessary to perform the customization once again. This is because we had to reset all toolbar settings for the changes to take place.

Below are pictures of [Chromatogram](#) and [Calibration](#) windows. Notice that the [Chromatogram](#) window has a dark yellow title bar, this is because it is opened from the second [Instrument](#). On the other hand, title bar of the [Calibration](#) window is green because it is opened from the third [Instrument](#). Windows opened from the fourth [Instrument](#) would have violet title bars.




Fig 7: Chromatogram window




Fig 8: Calibration window

Even outside of Clarity the icons for different file types have changed.




Fig 9: Windows OS - File Explorer - filetype icons

## 2.7 Batch Full version

- Batch dialog serves for reprocessing large number of chromatograms, calibration standards or measured sequence.
- We have added a field to select with which method you want to reprocess files - previous checkbox called *Reprocess by Instrument Method* was ambiguous as it was not clear which method was opened. In addition if you wanted to reprocess with a different method (than opened on the [Instrument](#)), you had to close [Batch](#) and open method on the [Instrument](#). Now simply choose file type, select reprocess and using the three dot button select by which method it will be reprocessed.
- Similar behavior applies to print, where it requires to select which *Report Style* will be used for the print.




Fig 10: Batch - Reprocess by Method and Print Results

## 2.8 Good Laboratory Practice (GLP)

GLP is for us the top priority and for this reason we are continually working on ways to ensure Clarity stays compliant ready. Based on experience from our customers and external auditors we have implemented the following features:

- Global Audit Trail
  - Located on a separate tab in the Station Audit Trail.
  - Logs critical stations operations: station starts/shut downs, creation of a new daily audit trail, changes in the configuration and user accounts, etc.
  - Logs are stored in one convenient file. Once it exceeds the allowed size, successive file is created.
- The option "*Disallow to Use Unsaved (Modified) Method*" has been removed from the [GLP Options](#). The reason is that change of method handling will not allow to *Send* and thus measure with a modified method.
- Improved security of password (\*.psw) file. This file is automatically created when setting up user accounts. Note that \*.psw file is not backward compatible, meaning that if it was created under Clarity 8.0, it cannot be used in older versions of Clarity.

## 2.9 MS Extension - enhancements ✓ Full version

Upon suggestions from our users, we have redesigned the **MS Search** which is an important part of the MS Extension. Different workflows were taken into account while developing the redesigned functionality.

### 2.9.1 Unified features

**Single Compound Search** and **Automatic Compound Search** possessed the same functionality but with different labels and location on the respective tabs. Shared functionality is now located within the group **Search Options** and under the same labels which should not cause any confusion.

The screenshot shows the MS Search dialog box with the following details:

- Search Options (highlighted in red):**
  - Min Match Factor: 500 (0 .. 1000)
  - Restrict m/z Range: From 40 To 250
  - Use Selected m/z
  - Search Only Selected
  - Search All But Selected
- MS Library Search Table:**

	Chrom. RT	Match	R. Match	Prob. [%]	Compound Name	Library	ID	Formula	MW	CAS No.
1	4,815	747	804	93,05	Methylene Chloride	MAINLIB	775	0 84	84	75-09-2
2	5,762	696	767	59,77	Pentane, 3-methyl-	MAINLIB	1301	0 86	86	96-14-0
3	6,845	749	749	86,84	Ethyl Acetate	MAINLIB	2029	0 88	88	141-78-6
4	7,168	692	755	63,64	Trichloromethane	MAINLIB	726	0 118	118	67-66-3
5	7,733	806	817	75,05	Benzene	MAINLIB	96	0 78	78	71-43-2
6	7,857	868	868	80,83	Pentane, 2,2,4-trimethyl-	MAINLIB	1579	0 114	114	540-84-1
- Chromatogram:** Shows Relative Intensity [%] vs m/z. A peak is labeled at 7,733 [min] (Spectral Data). A red line represents Benzene.

Fig 11: MS Search - Search Options group

### 2.9.2 Single Compound Search

**Single Compound Search** has undergone major enhancements. The **Match Spectrum at** label has been renamed to **Search In Ret. Time** ①. More importantly **MS Search** can now be performed over a range of retention times - simply check the **to** checkbox and fill the second retention time.

Now a specific m/z can be excluded from **MS Search** - check the **Use Selected m/z**, fill values which should be excluded from search and then select the **Search All But Selected** ②.

Support for subtraction background was implemented, in the *Background Subtraction* group (3), check one or two background checkboxes and fill the range. Note that if you fill a range that falls within the searched retention time you will be notified about it.

*Preview Spectrum in Library* (4) opens a [NIST MS Search 2.0](#) external utility to overview the set parameters - actual search is performed by the *Search* button.

Compound to be added to the MS Method is now selected using the checkbox on the specific row (5). Note that only ONE compound can be selected.

The screenshot shows the MS Search dialog box with the following details:

- Search In Ret. Time [min]:** from 6,397 to 7 (1)
- Search Options:**
  - Min Match Factor: 400 (0..1000), Max Hits: 3
  - Restrict m/z Range: From 35 To 150
  - Use Selected m/z: m/z 1..4: 129, 131, 164, 166
  - Search Only Selected
  - Search All But Selected (2)
- Background Subtraction (3):**
  - Background 1 [min]: from 5 to 7
  - Background 2 [min]: from to
  - Warning: Background overlaps with searched retention time!
- Buttons:** Search (4), Preview Spectrum in Library (4)
- MS Library Search Table:**

Match	R. Match	Prob. [%]	Compound Name	Library	ID	Formula	MW	CAS No.
1 <input checked="" type="checkbox"/>	914	914	98,09 Ethyl Acetate	Demo_ms	20	C <sub>4</sub> H <sub>8</sub> O <sub>2</sub>	88	141-78-6
2 <input type="checkbox"/>	702	702	98,09 Ethyl Acetate	MAINLIB	2029	C <sub>4</sub> H <sub>8</sub> O <sub>2</sub>	88	141-78-6
3 <input type="checkbox"/>	615	617	0,98 1,3-Propanediol, diacetate	MAINLIB	1689	C <sub>7</sub> H <sub>12</sub> O <sub>4</sub>	160	628-66-0
- Buttons:** Add Selected to Method
- Mass Spectrum Plot:** Shows Relative Intensity [%] vs m/z. Two series are plotted: 6,397 - 7,000 min (Spectral Data) in blue and Ethyl Acetate in red.
- Buttons:** Close, Help

Fig 12: MS Search - Single Compound Search

## 2.9.3 Automatic Compound Search

Apart from the changes found in the *Search Options* group, several enhancements have been implemented. The *Restrict Chrom. Retention Time* label has been renamed to *Search In Ret. Time* (6) and placed at the top of the tab - as it is one of the first steps that must be performed.

Search results from the MS library are still displayed in a table. However it is much easier to select multiple compounds to be added to the MS Method. Simply check the checkbox (7) on the row of the compound(s) that you want to add.

Selecting or deselecting all compounds can be now easily performed using the checkbox (8) below the table.

MS Search

Single Compound Search Automatic Compound Search Target Compound Search

Use Signal: TIC m/z: 35

Search In Ret. Time [min]:  Whole Chromatogram from 3 to 14,3

Search Options

Min Match Factor: 500 (0 .. 1000)

Restrict m/z Range From: 40 To: 250

Use Selected m/z m/z 1..4:

Search Only Selected

Search All But Selected

Search

Search in Libraries:

Demo.ms

MAINLIB

MS Library Search

	Chrom. RT	Match	R. Match	Prob. [%]	Compound Name	Library	ID	Formula	MW	CAS No.
1	4,815	747	804	93,05	Methylene Chloride	MAINLIB	775	0 84	0	75-09-2
2	5,762	696	767	59,77	Pentane, 3-methyl-	MAINLIB	1301	0 86	0	96-14-0
3	6,845	749	749	86,84	Ethyl Acetate	MAINLIB	2029	0 88	0	141-78-6
4	7,168	692	755	63,64	Trichloromethane	MAINLIB	726	0 118	0	67-66-3
5	7,733	806	817	75,05	Benzene	MAINLIB	96	0 78	0	71-43-2
6	7,857	868	868	80,83	Pentane, 2,2,4-trimethyl-	MAINLIB	1579	0 114	0	540-84-1

Select/Deselect All

Add All Selected to Method

Rel. Intensity [%]

7,733 [min] (Spectral Data)

Benzene

m/z

Close Help

Fig 13: MS Search - Automatic Compound Search

## 2.9.4 Independent size of views

Chromatogram graph and spectrum graph can now be resized independently, therefore you can customize it to your needs - move mouse cursor over the separation line and drag the line up or down, whether you want to display larger chromatogram graph or spectrum graph.




Fig 14: Chromatogram - split views



## 2.10 Internal Standards (ISTD) ✓ Full version

Upon suggestions, we have increased number of internal standards from 5 up to 10.

In order to make the most of the use of increased ISTD, it has been added to all relevant places:

- [Single Analysis](#)
- [Sequence](#)
- [Chromatogram](#)
- [Calibration](#)
- [Prints and Exports](#)




Fig 15: Single Analysis - Edit ISTD Amounts




Fig 16: Chromatogram - Edit ISTD Amounts

## 2.11 Calibration

### 2.11.1 Details of Calibration point

Now it is much easier to review what calibration points are being used for constructing the calibration curve.

It is accessible directly from the *Calibration Summary Table* - upon pressing the three dots button (on the row of the desired compound) in the column *Rec No.*, dialog *Details of Calibration* for the respective compound on the given level pops up. It is not necessary to go on the tab of the respective compound to view this information.




Fig 17: Calibration - Details of Calibration Point

## 2.12 Chromatogram

### 2.12.1 Enhanced error format

Enhanced error format for *Result Table* and *All Signals Result Table* mean that it is much easier to understand what is the problem and thus what must be repaired.

Each error/warning is located on a separate row, right above the tables.




Fig 18: Chromatogram window - error/warning above All Signals Result Table

### 2.12.2 Import Chromatogram

Command *File - Import Chromatogram...* opens *Import* dialog which now contains *Apply Method* field for method. Since method is an integral part of any chromatogram, removal of the method opened on the *Instrument* meant that we had to provide some other way for imported chromatograms to "have" method. For this reason the *Import* dialog has a new field where the user can choose a method that will be applied to the resulting chromatogram.

Parameters that will be copied from the method are: integration table, parameters from the calculation and advanced tabs. By default the

method that is set here is the same as selected in the [Single Analysis](#) dialog.




Fig 19: Import dialog - Apply Method

### 2.12.3 Mathematical Operations

Just as the *Import* dialog, so does the *Mathematical Operations* dialog contain a field for method. The functionality is the same as described above.




Fig 20: Mathematical Operations dialog - Apply Method

## 2.12.4 New icons for Set Model and Copy from Model

Functionality of *Set Model* and *Copy from Model* is commonly used and for this reason we have added icons which can be easily added (right mouse click - *Customize* - drag icon on the toolbar) on your toolbar, thus speeding up the work.




Fig 21: Chromatogram window - Method menu

## 2.13 Print

### 2.13.1 Report Setup - Word Wrap Long Texts

Compound names and other texts can be very long, therefore the option to *Word Wrap Long Texts* should be helpful in such situation. In the **Report Setup** on the *Results* tab, check the *Summary Table* checkbox and then check the *Word Wrap Long Texts*.




Fig 22: Report Setup dialog - Word Wrap Long Texts

The figure shows two side-by-side tables. The left table is a chromatogram table with columns: Sample ID, Sample, Sample Amount, and Reten. Time [min]. The right table is a report table with columns: Sample ID, Signal, Sample, Sample Amount, and Reten. Time [min]. Red boxes highlight the text truncation in the chromatogram and the full text in the report.

Sample ID	Sample	Sample Amount	Reten. Time [min]
Sample_Vial_6-1	Signal 1	Haloca	Sample 0,000 3,660
Sample_Vial_6-2	Signal 1	Haloca	Sample 0,000 3,667
Sample_Vial_7-1	Signal 1	Haloca	Sample 0,000 3,680
Sample_Vial_7-2	Signal 1	Haloca	Sample 0,000 3,673
Sample_Vial_8-1	Signal 1	Haloca	Sample 0,000 3,650
Sample_Vial_8-2	Signal 1	Haloca	Sample 0,000 3,643
Sample_Vial_9-1	Signal 1	Haloca	Sample 0,000 3,647
Sample_Vial_9-2	Signal 1	Haloca	Sample 0,000 3,660

Sample ID	Signal	Sample	Sample Amount	Reten. Time [min]
Sample_Vial_6-1	Signal 1	Haloca bars:	Sample 0,000	3,660
Sample_Vial_6-2	Signal 1	Haloca bars:	Sample 0,000	3,667
Sample_Vial_7-1	Signal 1	Haloca bars:	Sample 0,000	3,680
Sample_Vial_7-2	Signal 1	Haloca bars:	Sample 0,000	3,673
Sample_Vial_8-1	Signal 1	Haloca bars:	Sample 0,000	3,650
Sample_Vial_8-2	Signal 1	Haloca bars:	Sample 0,000	3,643
Sample_Vial_9-1	Signal 1	Haloca bars:	Sample 0,000	3,647
Sample_Vial_9-2	Signal 1	Haloca bars:	Sample 0,000	3,660

Fig 23: Cut off text in Chromatogram vs. full text in Report

### 2.13.2 Report Setup - Force Inverted

Sometimes printing of *Summary Table* is complicated as there can be more columns than rows. Clarity prints tables (for example from the *Results* tab) according to the columns as displayed on the screen. Option *Force Inverted* swaps columns for rows and vice versa. Thus enabling to view and print *Summary Table* with better readability.




Fig 24: Report Setup dialog - Force Inverted

### 2.13.3 Report Setup - new Lab. Header variable

New header variable "Instrument Name" is now available on the *Lab. Header* tab of the **Report Setup**. Inserting "%e" into the Lab. Header row, inserts name of the **Instrument** (as named in the **System Configuration**) into the report.




Fig 25: Report Setup dialog - variable Instrument Name

### 2.13.4 Renamed Report Style INSTRUMENT.STY

**Instrument** window used to have its own report style, as any other Clarity window. During the redesign we have removed the option to open **Report Setup** directly from the **Instrument** window. For this reason we have renamed the INSTRUMENT.STY to ANALYSIS.STY which is now the default print style for **Single Analysis** and **Sequence**.



## 2.14 Installation

### 2.14.1 Windows XP are no longer supported by Clarity 8.0

Microsoft has ended support for Window XP on April 8, 2014. For more than 4 years after the official end of support, we have allowed and supported Clarity on Windows XP. Now Clarity 8.0 can't be installed on any PC running Windows XP operating system.

### 2.14.2 Changes in the Typical installation

The **INT5**, **U-PAD**, **CB11** and **CB20** drivers (not supported on 64bit operating systems) were removed from the *Typical* Clarity Installation. In order to install them, select *Custom* and select required component(s) or alternatively select the *Full* installation.




Fig 26: Clarity Installation - Typical

## 2.15 Other changes

- [System Configuration - Available Control Modules](#) dialog includes a new filter "*Only Installed*" to display only modules that are installed on the PC.
- [LC Monitor](#) - for better understanding the buttons *Purge* and *Manual Flow* were renamed to *Manual Flow* and *Modify Gradient*.
- [Audit Trail](#) - better readability of the information displayed in the *Info* column.
- Command Line - message about restarting instrument will be suppressed.
- Clarity2Go - please update your Clarity2Go application to get a proper new instrument states which are now more self-explanatory.
- Teamviewer QS - tool for online support has been updated to version 13.
- DHA Extension - minor improvements in the DHA.CSV table.
- SST Extension - new option to restore default columns.
- Various known bugs have been fixed. See *What's new* in the [About](#) dialog of your Clarity.

## 3 New and updated control modules

This section contains new and updated control modules introduced in Clarity 8.0.

### 3.1 Advion

**Updated:**

- Advion Expression CMS - API updated to version 5.0.20.1.

### 3.2 CoruiTech

**New:**

- RainbowC HPLC - control module is now in the Testing state.
- Rainbow HPLC - control module is now in the Testing state.
- Rainbow detector - control module is now in the Testing state.

### 3.3 CTC

**Updated:**

- CTC PAL3 System driver updated to version 1.1.0.16.

### 3.4 Dani

**Updated:**

- Dani Master GC driver updated to version 1.6.9.0.

### 3.5 Ecom

**New:**

- Ecom ECF 2000 fraction collector - control module is now in the Testing state.

**Updated:**

- Flash12DAD and Flash06DAD - updated to version 2.4.0.0.
- ECD2000 - updated to version 3.3.0.0.
- ECDA2000 - updated to version 2.2.0.0.
- ECO2000 - updated to version 2.7.0.0.
- EcomECP2000 - updated to version 3.2.0.0.
- Panda30HID - updated to version 1.7.0.0.
- IOTA - updated to version 2.5.0.0.

### 3.6 HTA

**Updated:**

- HTA HT1500L autosampler driver updated to version 1.3.0.0.

### 3.7 KNAUER

**Updated:**

- Various KNAUER drivers updated to version 7.4.000.5067.

### 3.8 Shimadzu

**Updated:**

- Shimadzu AOC-20 sampler - added support for 0.5 and 5 µL syringe size.

### 3.9 Schambeck

**Updated:**

- Schambeck RI2000 detector - added support for models with firmware older than 4.86.

### 3.10 YL Instruments

**New:**

- ChroZen UHPLC ALS autosampler- control module is now in the Released state.
- ChroZen UHPLC CC thermostat - control module is now in the Released state.
- ChroZen UHPLC Pump pump - control module is now in the Testing state.
- ChroZen UHPLC UVD detector - control module is now in the Released state.

**Updated:**

- YL Instruments YL6500 GC driver updated to version 1.0.1.10.