

SW Protection - Clarity SW since version 2.4

Clarity Chromatography Software since version 2.4 (Clarity p/n C50, Clarity Lite p/n C40, Clarity Offline p/n C59) is protected by a dongle (hardware key, either USB or LPT). The earlier SW protection (CSW 17, CSW 32, and Clarity up to version 2.3) was based on a combination of a user code and a specific A/D converter.

Both distributors and customers benefit from higher flexibility. The essential difference and advantage of the dongle protection is the ability to use any A/D converter* with any Clarity SW. The simplified ordering procedure does not require specifying A/D cards for the purchase of Clarity. Distributors are able to keep stocks more easily now.

* Please refer to the compatibility table reflecting various operating systems.



Rockey LPT and USB hardware keys (dongles)

DataApex will distribute USB dongles by default. LPT dongles will be available upon request. Every dongle includes information on the S/N of Clarity and is accompanied by the corresponding user code.

IMPORTANT: Please note that a loss of a dongle (hardware key) is considered a loss of the license right to use Clarity SW.

A dongle installation is a part of the Clarity SW installation (this is done automatically). More details can be found in the *Getting Started* manual.

See the *Getting Started* manual (M006) at: <https://www.dataapex.com/downloads/26057/view>

Clarity ver. 2.4 release date is February 1st, 2005. All Clarity users may update to Clarity version 2.4 for a fee (refer to *Price list 2005*). Users of Clarity who ordered the license after August 1st, 2004, will receive this update free of charge (upon request).

For more information, contact sales@dataapex.com